

## 3-Wheel Cruiser

The GF 3-Wheel Cruiser is designed to provide convenient mobility for those individuals that do not require the extra stability of a 4-wheel rollator. This new light weight aluminum walking aid allows the user to have more freedom in movement to perform daily activities.

### Features

- Light weight aluminum frame
- Available in 3 new attractive colors: Pewter Metallic, Metallic Burgundy and Bondi Blue
- Ergonomic hand grips with locking loop brakes
- Adjustable handle height for broad range of users
- 8" black wheels for indoor and outdoor use
- Comes complete with zippered carry pouch
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0143
- Instructions included

<b>609201B</b>	Bondi Blue	1 ea/cs
<b>609201M</b>	Metallic Burgundy	1 ea/cs
<b>609201P</b>	Pewter Metallic	1 ea/cs



609201P

### Specifications:

<b>Overall Width:</b>	Rollator 28" Carry Pouch 14.17"
<b>Width Inside Hand Grips:</b>	14 ½"
<b>Width Inside Rear Legs:</b>	22"
<b>Height Adjustment:</b>	33" (min.) - 37" (max.)
<b>Patient Height:</b>	Approx. 5' 6" – 6' 5"
<b>Product Weight:</b>	13.8 lbs
<b>Weight Capacity:</b>	Rollator 250 lbs EVENLY DISTRIBUTED Carry Pouch 22 lbs
<b>Packaging:</b>	1 ea/carton
<b>Warranty:</b>	5 Year Limited Lifetime Warranty



609201B



609201M



GRAHAM-FIELD

[www.grahamfield.com](http://www.grahamfield.com)