

## Adjustable Height Versaframe

- Versaframe is designed to add support and enhance safety while toileting
- Contoured waterfall armrests are designed to add comfort and provide an excellent gripping surface
- Anodized aluminum frame is sturdy and lightweight
- Quickly secures to most toilets
- Arm height adjusts from 28" – 32"
- Arm width adjusts from 18" – 21"
- Easy tool-free assembly
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0243
- Lifetime limited warranty

<b>6460A</b>	Toilet Safety Rails, Versaframe	4/cs
<b>6460R</b>	Toilet Safety Rails, Versaframe, Unassembled, Retail Package	3/cs
<b>6465A-1</b>	Toilet Safety Rails, Versaframe, Unassembled, brown box	1 ea



### Specifications

Width Between Arms .....	18" - 21"
Overall Width .....	21.75"
Maximum Weight Capacity.....	250 lb
	EVENLY DISTRIBUTED
Replacement Tip.....	LX9141

