

Versa Height Raised Toilet Seat

- Features four, clip-on brackets that fit most toilets
- Angle and height are adjustable to maximize comfort and reduce bending and patient effort required
- Solid molded plastic seat ensures durability
- Versa Guard brackets will not scrape or mar toilets
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0244
- Lifetime limited warranty

6490A 1 ea
6492A with 2 rear locking brackets 1 ea



6490A



6492A



Specifications

Seat Width 14 ¾"
 Seat Depth 15 ½"
 Seat Height Range 3 ¾" - 5 ¾"
 Maximum Weight Capacity 300 lb EVENLY DISTRIBUTED