

Heavy Duty Bathtub Safety Rail

- Unique design provides a good gripping surface from any direction
- Heavy-gauge steel is powder-coated for rust resistance
- Provides two gripping surfaces to facilitate both entry and exit from tub
- Extra large suction cups add stability
- Fits tubs with walls 3" - 7 ¼" thick
- Not recommended for use with fiber glass tubs
- Color: White
- Attractive retail package
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Lifetime limited warranty

69603A 3/cs



GRAHAM-FIELD

www.grahamfield.com