

Tub-Guard® Bathtub Safety Rails

- Textured, contoured handle provides a good gripping surface
- Clamp holds firmly to tub wall without marring tub surface
- Lightweight, compact design is excellent for travel
- Color: White
- Fits tubs with walls 2 ½" - 6 ¼" thick
- Not recommended for use with fiber glass tubs
- Tool-free installation
- Available in two sizes: standard (12" height) and tall (16 ½" height)
- Tall model provides two gripping surfaces to facilitate both entry and exit from tub
- Attractive retail packaging
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Lifetime limited warranty

69623A Standard 12" 3/cs
69663A Tall 16 ½" 3/cs



69623A



69663A

