

Bariatric Steel Folding Commode

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame
- Welded steel construction provides strength and stability
- Folds flat for convenient storage and travel (to less than 7")
- Tool-free assembly
- Sturdy plastic seat snaps-on easily
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0168
- Maximum Weight Capacity 650 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty

7109A-2 2/cs



Specifications

Width Between Arms	24 ½"
Overall Width	30"
Seat Width	14"
Seat Depth	16 ½"
Seat Height Range	17 ½" - 21 ½"
Maximum Weight Capacity	650 lb EVENLY DISTRIBUTED

