

Shower Transport Chair

- Designed for safe transport to and from the shower
- Anodized, non-corrosive aluminum frame construction will not rust
- Comfortable nylon backrest removes for easy cleaning
- Cushioned armrests enhance user comfort
- 5" swivel casters with rear wheel safety locks facilitate maneuverability in tight spaces
- Rolls easily over most toilets
- Allows complete access to body for easy cleaning during bathing
- Maximum Weight Capacity 250 lb EVENLY DISTRIBUTED
- · Limited lifetime warranty

7910A-1 1 ea

Specifications

Width Between Arms	17 ½"
Overall Width	21 ½"
Overall Depth	21 ½"
Overall Height	38"
Seat Width	14"
Seat Depth	16 ½"
Seat Height	22"
Maximum Weight Capacity250 lb EVEN	NLY DISTRIBUTED



