



Manual Hi-Lo Flex Bed Shown with Spring Mattress Deck, Stabilizer Legs, and Stock Medium Oak Head Board and Foot Board.

Flex Bed™

- Height travel range: 14.5" to 26"
- Safe work load: 450 pounds*
- Hi-Lo Manual
- Spring sleep deck
- Stabilizing legs standard
- Durable powder-coated epoxy finish
- UL 2601-1 and CAN/CSA C22.2 No. 601-1-M90 approved
- Warranty:
 - 15 years on frame
 - 2 years on bed mounted accessories
 - 1 year Head/foot boards

* including resident, support surface, head/foot boards, and accessories

76" LENGTH	80" LENGTH	Mattress Deck Surface	Bed Type	Units
AF4000	AF4100	Spring Deck	Manual Hi-Lo	1 ea

* Bed, head and foot boards, and rails sold separately.

Wall Saver Bumpers

A65476	Wall Saver Bumpers, 76" Bed
A65480	Wall Saver Bumpers, 80" Bed

Options and Accessories

Factory-Installed Options

999-0700-353	Pull Pin Lock Out Pendant
554-0658-000	Edema Foot Ratchet
554-0769-000	Low Profile Matt Retainer Wire Form Kit

Head and Foot Boards*

AB53300	Rectangular, Solar Oak Finish
AB48300	Rectangular, Medium Oak Finish
AB80300	Rectangular, Montana Walnut Finish
AB29300CCR	Center Crown, Windsor Mahogany Finish
AB04300CCR	Center Crown, Wild Cherry Finish

* All boards above have T-mold edges.
Other styles and finishes available via custom order.

Side Rails and Assist Bars

A8340	Half Head Rails (Set)
A7050	Half Foot Rails (Set)
A7070GREY	Full Length Rail- 72" Length (Set)
AZ7350GREY	Assist Bar (1 Each)

Accessories

555-0206-902	Emergency Crank**
554-0639-000	Trapeze Adapter
556-0005-000	Traction Frame Adapter (2 Req)
554-0143-000	Trapeze Bar (Patient Helper Bar)
554-0190-000	Panel-Mounted I.V. Holder

800.554.9215 Fax 800.524.4113
View these and more products at www.grahamfield.com