

First Aid Kit - 50 person, Plastic Case

- Designed for small businesses and work site
- Suitable for up to 50 employees
- Includes assortment of bandages, gauze pads, tape, ointment and more
- Water-resistant plastic case with handle and can be wall mounted
- *Medicine & Antiseptics*
 - 30 BZK Antiseptic towelettes
 - 10 First aid/burn cream packs, 0.9 g, 1 box, 10/bx
 - 10 Antibiotic ointment pack, .9gm, 1 box, 10/bx
 - 4 Ammonia inhalant ampoules
- *Bandages*
 - 125 1" x 3" Adhesive plastic bandages
 - 10 2" x 4" Elbow/Knee bandages
 - 6 Knuckle fabric bandages
 - 6 Fingertip fabric bandages
- *Injury Treatment*
 - 3 5" x 9" Trauma pad
 - 8 3" x 3" Gauze dressing pads, 4 2-pks
 - 1 40" x 40" x 56" Triangular sling/bandage, w/2 safety pins
 - 1 1/2" x 10 yd. First aid tape roll
 - 4 Sterile eye pads
 - 2 3" Conforming gauze roll bandage
 - 3 2" Conforming gauze roll bandage
 - 8 2" x 3" Non-stick pads
 - 1 1 oz eyewash
 - 2 Tourniquets, latex free
 - 1 Finger splint
 - 1 4" x 5" Instant cold compress
 - 1 6" x 9" Instant cold compress
- *Reference / Instruments*
 - 1 CPR mask, 2 exam quality gloves, 3 BZK antiseptic towelettes, 1 set/bx
 - 1 First Aid Guide
 - 4 Exam quality vinyl gloves
 - 1 4-1/2" Scissors, nickel plated
 - 1 4" Tweezers, plastic, one time use
 - 1 Plastic case
- **Total 253 pieces**

**Meets OSHA recommendations
and ANSI/ISEA Z308.1-2009**

1799-50P Plastic case 10¾" x 11¼" x 3" 1 ea



GRAHAM-FIELD

www.grahamfield.com

© 2012 GF Health Products, Inc.
Atlanta, Georgia 30360 GF1200114RevA12