

Reachers

Lightweight aluminum reaching aids help avoid bending down and climbing on objects. A squeeze of the trigger closes jaws for a tight grip. Plus a magnetized tip helps to retrieve small metallic items.

Ideal for:


- Reaching in tight spaces such as behind furniture
- Grabbing object just out of reach in high spaces
- Picking up objects off the ground

5689	Basic Reacher	1 ea
5690	Standard Reacher	1 ea
5691	Deluxe Reacher	1 ea

	Basic Reacher (Item# 5689)	Standard Reacher (Item# 5690)	Deluxe Reacher (Item# 5691)
Helps avoid bending and climbing	✓	✓	✓
Magnetized tip for retrieving small metallic items	✓	✓	✓
Attractive retail package	✓	✓	✓
Extends Reach	32"	32"	30"
Jaw Width	3"	3"	2 ¾"
Product Lightweight	½ lb	½ lb	½ lb
Lifting Weight	2 lb	2 lb	2 lb
Folds for easy storage and transport	✓	✓	
Jaw mechanism	Saw tooth surface	Suction cups	Saw tooth & anti-slip surface
Hook to aid with getting dressed			✓
Jaws can be turned 360° to avoid awkward wrist turns			✓




Magnetized Tip



Deluxe version has 360° rotation




Basic and Standard versions fold for easy storage and transport



Basic



Deluxe



Standard



GRAHAM-FIELD

www.grahamfield.com