

## Wallmax™ Aneroid Sphygmomanometer

- Luminescent dial
- Large, 6", easy-to-read gauge
- Wall mount bracket with 110° swivel range
- Accumax™ nylon cuff with 8' coiled tubing
- Cuff can be stored behind gauge
- One (1) year warranty for manufacturing defects only

**V222B** Adult, **Latex Free** bladder Black Cuff 1 ea

**CR222** Crystal Faceplate 1 ea

**222** Adult, **Latex** bladder Blue Cuff 1 ea

**Caution: This product contains natural rubber latex, which may cause allergic reactions.**

Cuff Class	Limb Range (Inches)	Bladder (Inches)	Cuff Dimensions (Inches)*
Newborn	3.75" - 5.25"	1" x 3.55"	7" W x 2" H
Infant	5" - 7.5"	2.35" x 4.72"	10.5" W x 3" H
Child	7" - 10.5"	3.55" x 5.90"	13.5" W x 4" H
Adult	9.5" - 15"	4.72" x 8.66"	20.5" W x 5.5" H
LG Adult	13" - 19"	5.71" x 12.41"	25" W x 6.5" H
Thigh	17" - 26"	7.09" x 14.57"	29.5" W x 8.5" H

\* These are subjective measurements and can vary  $\pm 0.25"$



CR222



222



V222B

