

Fashion Folding Canes

- Available in a variety of attractive finishes
- Cane folds to an 11 inch length for convenient transport and storage
- Derby-style handle is contoured for comfort
- Easy push-button height adjustment with locking collet for added security
- Height adjusts from 33" – 37" in 1" increments
- 200 lb maximum weight capacity
- Designed to meet the requirements of HCPCS code: E0100
- Limited lifetime warranty

5950DF Combo Pack 12/cs
2 Golden Leopard
2 Paisley Purple
2 Paisley Green
3 Paisley Blue
1 Pink Flower
1 Purple Flower
1 Yellow Flower



GRAHAM-FIELD

www.grahamfield.com

© 2013 GF Health Products, Inc.
Atlanta, Georgia 30360 GF0700053RevC13