

Fashion Offset Canes

- Available in a variety of attractive finishes
- Offset handle places the user's weight directly over the base for maximum balance and control
- Easy push-button height adjustment with locking collet for added security
- Height adjustment 31"- 40"
- 250 lb maximum weight capacity
- Designed to meet the requirements of HCPCS code: E0100
- Limited lifetime warranty
- Soft, warm and contoured Nitrile grip for maximum comfort
- Sculpted grip reduces the chance of slippage, even when wet
- Comes complete with convenient wrist strap

5950B **Combo Pack** 6/cs
3 Paisley Green
3 Paisley Blue

