


## Palm Aneroid Sphygmomanometer

- Ambidextrous rotating bulb and valve
- Nylon cuff with index & range markings
- Latex-free one (1) tube bladder
- Durable aluminum manometer
- Large and easy-to-read 65 mm dial
- Ribbed inflation bulb for easy and secure grip
- Color coordinated gauge and cuff, grey or orange
- Deluxe push-button air release valve
- Deluxe end valve reduces dust build-up
- Twenty (20) year calibration warranty
- Deluxe nylon carry case is included



- 215 Adult, Grey 1 ea  
 215C Child, Grey 1 ea  
 215X Large Adult, Grey 1 ea  
 216 Adult, Orange 1 ea



<u>Cuff Class</u>	<u>Limb Range (inches)</u>	<u>Bladder (inches)</u>	<u>Cuff Dimensions (inches)*</u>
Newborn	3.75" - 5.25"	1" x 3.55"	7" W x 2" H
Infant	5" - 7.5"	2.35" x 4.72"	10.5" W x 3" H
Child	7" - 10.5"	3.55" x 5.90"	13.5" W x 4" H
Adult	9.5" - 15"	4.72" x 8.66"	20.5" W x 5.5" H
LG Adult	13" - 19"	5.71" x 12.41"	25" W x 6.5" H
Thigh	17" - 26"	7.09" x 14.57"	29.5" W x 8.5" H

\* These are subjective measurements and can vary  $\pm 0.25"$