

## Multi-Cuff-5 Blood Pressure Kit

- Ambidextrous rotating bulb and valve
- Orange nylon latex-free cuffs with index & range markings
- Five (5) different sizes: infant, child, adult, large adult and thigh
- Durable aluminum manometer
- Large and easy-to-read 65 mm dial
- Latex-free inflation one (1) tube bladder
- Ribbed latex-free inflation bulb for easy and secure grip
- Deluxe push-button air release valve
- Deluxe end valve reduces dust build-up
- Twenty (20) year calibration warranty
- Padded carry case with extra storage compartment - hook and loop closure

206 Orange 1 ea



Cuff Class	Limb Range (inches)	Bladder (inches)	Cuff Dimensions (inches)*
Newborn	3.75" - 5.25"	1" x 3.55"	7" W x 2" H
Infant	5" - 7.5"	2.35" x 4.72"	10.5" W x 3" H
Child	7" - 10.5"	3.55" x 5.90"	13.5" W x 4" H
Adult	9.5" - 15"	4.72" x 8.66"	20.5" W x 5.5" H
LG Adult	13" - 19"	5.71" x 12.41"	25" W x 6.5" H
Thigh	17" - 26"	7.09" x 14.57"	29.5" W x 8.5" H

\* These are subjective measurements and can vary  $\pm 0.25"$