

Walkabout Junior Four-Wheel Rollator

This light-weight Walkabout Junior rollator is ideal for the pediatric user, has the same great features as the Walkabout Lite, weighs only 11.5 lbs, and is available in two attractive colors. It folds quickly and easily into a compact unit for storage or transport. Its features include a padded seat, a backrest for extra comfort, ergonomic hand grips, and easy-to-operate hand brakes.

- Epoxy-coated aluminum frame
- Locking loop brakes
- Removable wire storage basket
- Adjustable handle height
- 6" wheels for indoor or outdoor use
- Designed to meet the requirements of HCPCS codes: E0143 + E0156
- Includes laminated color operating instructions
- Limited lifetime warranty

RJ4301R Burgundy 1 ea

Specifications

Width Inside Hand Grips.....	12"
Overall Width.....	20"
Overall Depth.....	24"
Seat Width.....	12"
Seat Depth.....	12"
Seat Height.....	18.5"
Caster Diameter.....	6"
Minimum Handle Height.....	27.5"
Maximum Handle Height.....	32"
Maximum Weight Capacity.....	300 lb

EVENLY DISTRIBUTED



Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.