

## Lab Chair - High Bench

### Advantages:

- Durable Polyurethane
- Adjustable height and tilt backrest
- Height range: 25.5" – 35.5"

### Features:

- Polyurethane seat and backrest
- Chrome-plated foot ring
- Chrome-plated gas cylinder
- Chrome-plated base
- Nylon casters

**1932LH** 1 ea

### Specifications:

Weight Capacity .....300 lb - EVENLY DISTRIBUTED  
Height Range.....25.5" – 35.5"  
Seat / Backrest ..... 16.9" W x 16.5" D / 15.35" W x 11"H  
Base ..... 24.4" Ø  
Caster ..... 2.5" Ø

### Shipping Info:

Gross Wt.....40.35 lb  
Net Wt.....35.94 lb  
Carton.....26.6" W x 25.2" L x 8.7" H

