

## PVC Shower Commode Chair

- Constructed of the highest quality grade PVC
- Deluxe elongated heavy-duty open front seat
- Highest quality casters are made to last and roll easily over surfaces
- **89180** comes standard with 3" twin wheel casters; 2 locking and 2 non-locking
- **89220** comes standard with 3" single wheel casters; 2 locking and 2 non-locking
- Quick drying, removable back support is fluid resistant, anti-bacterial and anti-fungal
- One year limited warranty on PVC frame; six month limited warranty on all other materials



*Highest quality casters are made to last*



*Contoured flip-up seat for easy cleaning*

**89180** 18" 1 ea  
**89220** 22" 1 ea

<b>Specifications</b>	<b>89180</b>	<b>89220</b>
Overall Height.....	41.3"	41.3"
Seat Height.....	19.5"	19.5"
Overall Width.....	21.5"	25.5"
Overall Depth.....	21.5"	21.5"
Seat Depth.....	18.3"	18.3"
Width Between Arms.....	18"	22"
Weight Capacity.....	300 lbs	375 lbs
	(EVENLY DISTRIBUTED)	



*Comfortable fluid resistant, anti-bacterial and anti-fungal removable back support*

