

**Medicopaste® Bandage**

- Application – Aids in compression treatment:
  - Lymphatic edema
  - Thrombophlebitis of the superficial saphenous vessels
  - Eczema of the leg
  - Ulcer of leg
- Can be used for certain static conditions of the leg and foot, such as post-fracture foot conditions, spastic weak foot conditions
- Can function as a “flexible cast” in treatment of sprains, strains, and dislocations
- Contraindication in arterial vascular diseases
- Zinc gelatin premium gauze bandage
- Made in USA, GH Health Products, Inc Rhode Island Facility
- ISO 13485:2003 Certified Facility

**1565 3**    3" x 10 yds    1/dz  
**1565 4**    4" x 10 yds    1/dz

*With Calamine*

**1565C-3**    3" x 10 yds    1/dz  
**1565C-4**    4" x 10 yds    1/dz

