

## Grafco® Deluxe Exam Lamps Series

- **Maximum of 60 Watt** or 13 Watt CFL equivalent
- 115/120V AC. UL listed 3-wire plug with 8ft grounded cord
- Flared, doubled-wall, dual cool ventilated reflector
- Heavy gauge anodized aluminum gooseneck shade
- One piece phenolic socket with back-turn knob switch
- Height range: 34 ½" – 53½"

*Exam Lamps with 3 wire 8 ft. grounded cord/plug.*

<b>1698</b>	Chrome-plated Base, Aluminum Shade	1 ea
<b>1698C</b>	With Clutch-collar Lock	1 ea
<b>1698CM</b>	With Mobile Base/Clutch-collar Lock	1 ea
<b>1698M</b>	With Mobile Base	1 ea

*Exam Lamps with 3 wire 8 ft. grounded cord/hospital grade plug.*

<b>1698-1</b>	Chrome-plated Base, Aluminum Shade	1 ea
<b>1698-1C</b>	Chrome-plated Base with Clutch-collar Lock	1 ea
<b>1698-1CM</b>	With Mobile Base/Clutch-collar Lock	1 ea

