

# Gowllands Ophthalmoscope/Otoscope Set

## Gowllands Ophthalmoscope/Otoscope Set

- Requires 2 AA batteries (not included)
- Sterilizable – autoclavable
- Made in United Kingdom

### *Ophthalmoscope:*

Length: 7.75"

Power lens range: from +25D to -40D

### *Otoscope:*

Length: 6.75"

Four aural specula

One nasal speculum power of 3.5

Magnifying lens –

**5100X** 1 ea



GRAHAM-FIELD

[www.grahamfield.com](http://www.grahamfield.com)

© 2012 GF Health Products, Inc.  
Atlanta, Georgia 30360 GF1200077RevA12