

Folding Set Up Screens

- 6ml white vinyl panels are flame-retardant*
- Easy to set up and to dismantle for storage
- Vinyl panel is washable and replaceable on snap-out rods
- Frame is 1" (25mm) anodized aluminum tubing
- 2" (51mm) twin wheel castors allow for easy movement

*NFPA 701 Test Method 1-2010 Edition

4297W	Three section	1 ea
4299C	Replacement Casters	1 set
4299W	Replacement White Panel	1 ea

Specifications

Each Panel: 27.5" x 56" (700 mm x 1,420 mm)

Overall Height: 70.8" (1,800 mm)

Weight: 6.75 kg

