

Funnels

- Made of linear polyethylene
- Stackable and dishwasher safe
- Sterilization Guidelines

Boiling:

100°C or 212°F, 15 minutes maximum

Gas: Ethylene Oxide

30°C – 60°C

86°F – 140°F

There should be no distortion with Ethylene Oxide Gas Sterilization temperature range.

Please allow ample time for the funnels to dissipate the absorbed gas after sterilization.

Note: Sterilization guidelines were provided by the raw material supplier, Dow Chemical Company.

3618A	2 oz., Top size: 2 ¾"	1 ea
3618C	8 oz., Top size: 4 ¼"	1 ea
3618D	16 oz.	1 ea
3618E	32 oz.	1 ea

