

Physician Mechanical Beam Scale with Wheels

- Dual measurement in LB and Kg
- Weight capacity up to 490 lb (220Kg), using counterpoise weight of 50 lb (20 Kg)
- Weight can be read from the front or back
- Retractable aluminum height rod, Range of 24 inch – 84 inch (60 cm – 210 cm)
- Heavy-duty roller wheels allow easy portability of the scale
- Durable enameled steel body
- Cast Aluminum beam
- Removable slip-resistant platform cover
- No parts contain Latex: Latex Free
- **Caution:** Placing the scale on a soft surface such as carpet is not recommended. This could result in inaccurate weight measurement. Ensure scale is on a level, stable and flat surface when in use.
- Meeting the OIML, International Organization of Legal Metrology, International Recommendation R76

Warranty:

Two (2) year limited warranty on manufacturer defects. This warranty does not include calibration of the scale. See User Instructions for the details.

HT485 1 ea

Specifications

Max. Weight Capacity:	490 lb x 0.25 lb (220 kg x 0.1 kg)
Height Rod Range:	24" – 84" x 0.25" (60 cm – 210 cm x 5 mm)
Platform Dimensions:	14.76" x 10.83" (375 mm x 275 mm)
Overall Dimensions (L x W x H):	20.87" x 10.83" x 58.46" (530 mm x 275 mm x 1,485 mm)
Net Weight:	35.5 lb (16 Kg)
Gross Weight:	45.2 lb (20.5 Kg)
Cubic Dimensions:	59.5" x 20.7" x 6.3" (151 cm x 52.5 cm x 16 cm)

Specifications are subject to change without notice.

