

### Features:

- 3" threaded stem, heavy duty casters
- Floor plate to support liner
- Zippered mesh liner available in 3 anti-bacterial mesh or anti-bacterial solid vinyl colors
- 28 gallon capacity per bag
- Push handle for easy transport

### Optional leak-proof liner

Available in light blue, pink, mint green, and white  
 Super washable, sonically welded, meets OSHA regulations and NFPA701 flame retardant standards. Rot and mildew resistant.  
 Elastic flip-top.



Item	Description	Height	Width	Depth
8510-19	Single Bag	38 1/2"	18 3/4"	23 3/4"
8511-29	Double Bag	38 1/2"	18 3/4"	34 1/4"
8512-19	Single Bag, with foot pedal	38 1/2"	22 1/2"	23 3/4"
8517-29	Double Bag, with foot pedal	38 1/2"	22 1/2"	34 3/4"
8518	Optional leak-proof liner			

### Standard Mesh and Solid Vinyl Colors



Royal Blue



Forest Green



Maive