

- Provides transport, transfer, and commode functions all in one unit
- Arms drop independently below seat level for lateral transfers
- Cushioned armrests enhance user comfort
- Swing-away footrests can be extended or detached
- 5" rear swivel, locking casters facilitate maneuverability in tight spaces
- Padded seat cover
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Comes unassembled
- Designed to meet the requirements of HCPCS code: E0165
- Limited lifetime warranty



- 6810A** Includes 7 qt commode pail 1 ea
6810B 4 locking casters 1 ea

Specifications

Width Between Arms.....	18 ¼"
Overall Width	18 ½"
Padded Seat Width	16"
Padded Seat Depth	16"
Commode Seat Width	17 ¾"
Commode Seat Depth.....	17 ½"
Commode Seat Height.....	19 ½"
Maximum Weight Capacity	300 lb EVENLY DISTRIBUTED