



Traveler® HD 3G01 Series

The Traveler HD comes with a reinforced frame, heavy-duty 5/8" diameter axles, and a long lasting, embossed back and seat upholstery with a double inner liner for durability designed to support up to 500 lb users - EVENLY DISTRIBUTED. Available in 20", 22" and 24" seat widths, with 18" seat depth (1" longer than standard). Meets Medicare K0007 specifications.

Standard Features

- Hammer tone coating provides a long lasting and durable finish
- Chart pocket standard on back upholstery
- Industry compatible front rigging
- Durable aluminum black footplates

Warranty

5 year limited warranty



Dual Crossbrace

Adjustable Axle

Custom colored upholstery is the ideal solution to assist with inventory control, floor identification and theft prevention.



Black Gray Orange Bur gundy



Rehab Shower Commode 1202 Series

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. This chair features an overall narrow width and a corrosion-resistant stainless steel frame. The chair can be positioned over a standard toilet or be used as a portable, self-contained commode. Has a maximum weight capacity 250 lb - EVENLY DISTRIBUTED. Meets Medicare Code E0240 specifications.

Standard Features

- Completely removable full length arms for ease of transfer
- Recessed seat and frame side contours allow easy access for toileting activities while seated
- Removable, hygienic plastic pail and lid are included
- Non-corrosive composite footplates

Warranty

2 year limited warranty



Padded Seat, Oval Cutout



Padded Seat, Open Front (w/Hooks) -Optional

Visit www.grahamfield.com for a complete line of wheelchair accessories and loss prevention options.



GRAHAM-FIELD

GF Health Products, Inc. Atlanta, Georgia 30360
770.368.4700

Graham-Field is not responsible for typographical errors. Packaging, warranties and products are subject to change without notice. © 2014 Everest & Jennings, is a registered trademark of GFHealth Products, Inc.

GF1200138RevC15



Radiology

Orthopedic



Traveler® HTC 3J01 Series

The Traveler HTC (Hospital Transport Chair) provides the ideal solution to the hospital/clinical transport need by incorporating special attendant and inventory control features with a heavy duty frame reinforced at all critical load-bearing points. Available in 20", 22" and 24" seat widths, 18" seat depth (1" longer than standard). Meets Medicare K0007 specifications.

Standard Features

- Includes a fixed anti-theft post ideal in hospital/clinical settings
- Embossed black leatherette upholstery with double inner lining for durability
- Hammer tone coating provides a long lasting and durable finish
- Anti-tippers standard
- Chart pocket standard on back upholstery
- 24" x 1" black plastic wheels with black handrims and non-marking polyurethane tires
- Durable aluminum black footplates
- Maximum weight capacity 500 lb - EVENLY DISTRIBUTED

Warranty


5 year limited warranty



Push Handle/Anti-Fold Bar



Rear Mounted Attendant Brake



Anti-Theft Nuts



GRAHAM-FIELD



Traveler® SE 3E01 Series

The Traveler SE is the standard wheelchair of choice. Designed with dual rear axles and adjustable casters, this wheelchair can convert from a K0001 Medicare code to a K0002 Medicare code. The Traveler SE has industry compatible front rigging and has vinyl upholstery constructed with a strong inner lining for added support.

Standard Features


- Hammer tone coating provides a long lasting and durable finish
- Two position axle and caster fork provides 2" seat height adjustment
- 250 lb weight capacity - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002

Warranty


1 year limited warranty



Chart Pocket



Composite Footplate



Custom Embroidery Available
(minimum quantities apply)



Advantage® Recliner Wheelchair 3K01 Series

The Advantage Recliner comes with a higher back and removable full support headrest. The headrest includes an additional neck support that adjusts easily to fit every need. Anti-tippers and an anti-fold bar come standard. Meets Medicare K0001 / E1226 (Recliner Back) E0966 (Head Rest Extension) specifications for wheelchairs.

Standard Features

- Detachable desk arms
- Hammer tone finish provides durability and style
- Low maintenance black vinyl upholstery
- Available in 16", 18" & 20" seat widths, 17" seat depth.
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED

Warranty

5 year limited warranty



Head Rest



Anti-Fold Bar



Reclined Position



Traveler® LX 3E020 Series

The Traveler LX incorporates the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction and industry standard front rigging. The premium vinyl upholstery is reinforced with an extra heavy duty inner liner. This wheelchair comes in a chrome finish and provides the most effective combination K0001/ K0002 wheelchair available.

Standard Features


- Blue premium vinyl upholstery
- Chart pocket on back
- "Tool Free" adjustable front rigging
- Chrome finish provides a long lasting and easy to clean finish
- Industry compatible front rigging
- Gray non-marking tires
- Two position axle and caster fork provide 2" seat height and adjustment
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0001 and K0002

Optional Features

- Rear Mount Anti-tippers

Warranty

2 year limited warranty



Premium Vinyl Upholstery

Gray non-marking tires

3/4" Premium Polyurethane Foam

Extra Duty No Stretch Liner



Traveler® L4 3F02 Series, 3F03 Series for 18" seat depth models

The Traveler L4 is a high-strength lightweight K0004 Medicare coded wheelchair. The Traveler L4 has dual axles and multi-position caster forks for an adjustable seat height of 18 or 20 inches. The wheelchair comes standard with non-marking polyurethane tires, anti-tippers and a chart pocket.

Standard Features

- Weighs less than 34 lbs. without front rigging
- Black nylon upholstery with height adjustable back
- Industry compatible front rigging
- Flip back arms for lateral transfer
- Two position axle and caster fork provides 2" seat height adjustment
- 250 lb weight capacity (300 lb on models with 18" seat depth) - EVENLY DISTRIBUTED

Warranty

5 year limited warranty



Adjustable Armrest option



Adjustable Front Rigging



Compatible Front Rigging



Flip Back Arms