

- Sturdy 1" aluminum tubing provides maximum strength while remaining lightweight
- Dual-release folding mechanism enables user to fold walker sides independently
- 5" front fixed wheels
- Complete with glide tips that replace rear leg rubber tips to enhance walker maneuverability on all surfaces
- Designed to meet the requirements of HCPCS code: E0143
- One (1) year limited warranty



Specifications

Width Inside Handgrips.....	17"
Width Inside Rear Legs	20"
Overall Width.....	25 1/2"
Depth Front to Rear.....	20 1/4"
Depth Folded.....	4"
Minimum Height	32"
Maximum Height	39"
Weight Capacity	300 lb

*Measurements taken with walker in lowest position.
Specifications subject to change without notice.*

716270B-2	Adult	Blue	2/cs
716270BK-2	Adult	Black	2/cs
716270G-2	Adult	Gold	2/cs
716270P-2	Adult	Plum	2/cs
716270PK-2	Adult	Pink	2/cs
716270R-2	Adult	Red	2/cs



Color Options