

Features:

- Combination of folding walker and rising aid all in one
- Secondary handles provide stable assistance from a seated to a standing position
- Attractive, upscale Onyx finish
- Single release folding mechanism is designed to aid users with limited hand dexterity
- Can be used as a portable toilet safety frame
- Durable aluminum tubing provides strength while remaining lightweight
- Adjustable height fits a broad range of user heights
- 400 lb weight capacity
- Contoured, plastic grips for enhanced comfort and a secure hold
- Designed to meet the requirements of HCPCS code: E0135
- Includes laminated color operating instructions
- Three (3) year limited warranty



Specifications:

Width Inside Handgrips.....	18"
Width Inside Rear Legs.....	23"
Overall Width.....	26"
Depth Front to Rear.....	20"
Depth Folded.....	4.5"
Minimum Height.....	32"
Maximum Height.....	39"
Approx. User Height.....	5'4" - 6'2"
Product Weight.....	5.5 lbs
Weight Capacity.....	400 lbs

EVENLY DISTRIBUTED

700175C-2 2/cs
 700175CR Unassembled, Retail Packaging 1 ea



700175CR Unassembled version