

Features:

- Increased width and depth accommodates larger users
- 600 lb maximum weight capacity - EVENLY DISTRIBUTED
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0148
- Three (3) year limited warranty

Specifications:

| | |
|-------------------------|---------------|
| Width Inside Handgrips | 20" |
| Width Inside Rear Legs* | 23" |
| Overall Width* | 26 1/8" |
| Depth Front to Rear* | 18" |
| Height Range | 31" - 38" |
| Approx. User Height | 5' 4" - 6' 2" |
| Weight Capacity | 600 lb |

* Measured at maximum height



604070A 2/cs

with 5" wheels

604070W 2/cs