

- Single release folding mechanism is designed to aid users with limited hand dexterity
- Durable 1" aluminum tubing provides maximum strength while remaining lightweight
- Raised H-frame design facilitates over-toilet positioning
- Adjustable height fits a broad range of user heights
- Accommodates most Lumex accessories and leg attachments
- Designed to meet the requirements of HCPCS code: E0135
- Includes laminated color operating instructions
- Handgrip is contoured and textured to provide a better gripping surface
- Three (3) year limited warranty



*Handgrip is contoured and textured to provide a better gripping surface.*

### Specifications:

Width Inside Handgrips .....	17 1/4"
Width Inside Rear Legs* .....	19 3/4"
Overall Width* .....	23 1/4"
Depth Front to Rear* .....	19 1/2"
Depth Folded .....	4 3/8"
Height Range .....	32 1/2" – 39 1/2"
Approx. User Height .....	5' 4" – 6' 2"
Weight Capacity .....	300 lb
	EVENLY DISTRIBUTED



**613070A**    Adult    4/cs