

Features:

- Non-folding design provides extra rigidity for use where space is not an issue
- Angled legs provide maximum patient support in all directions and help resist tipping
- Angled side brace helps support patients as they stand up and sit down
- Fits easily around standard toilet fixtures
- Durable 1" aluminum tubing provides maximum strength while remaining lightweight
- Adjustable to fit a broad range of user heights
- Accommodates all leg attachments
- Designed to meet the requirements of HCPCS code: E0130
- Three (3) year limited warranty

Specifications:

Width Inside Hand Grips	17"
Width Inside Rear Legs	19 3/4"
Overall Width	22 1/2"
Depth Front to Rear*	20 1/2"
Minimum Height	31"
Maximum Height	38"
Approx. User Height	5' 3" - 6' 2"
Weight Capacity	400 lb

EVENLY DISTRIBUTED



1925A 4/cs