

Features:

- Encourages independence and self esteem by allowing freedom of movement and exercise
- Restraint free opening/closing gate with a safety knob
- Mesh available in 13 colors
- Adjustable safety belt to help prevent injuries
- Solid base seat with white cushion for extra comfort
- Constructed of the highest quality health care grade PVC
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED

Royal Blue mesh will be provided for online orders.

If you desire another color, please order via Customer Service at 1-800-347-5678.


Standard Mesh and Solid Vinyl Colors



Royal Blue



Forest Green



Mauve

Specifications:

Internal Width	18"
External Width	30.5"
Seat Height.....	20.25" - 24.75"
Height Adjustable Arms.....	28.75" to 33.25"
Total Height	40.75"
Total Length.....	35"
Walking Space.....	15" (Front of Seat to Gate)
Maximum Weight Capacity.....	300 lbs
	EVENLY DISTRIBUTED



80500

1 ea