

Teaching Stethoscope

- Excellent acoustical clarity.
- 22" latex-free tubing with connector allows both student and instructor to listen simultaneously.
- Anodized aluminum combination dual head chestpiece and chrome-plated brass rotating stem.
- Adjustable chrome-plated external spring binaural.

540

Black

1 ea



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com