

## Features:

- Provides assistance to those who have difficulty stepping over bathtub walls
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Drainage holes in seat minimize water build-up
- Side support arm aids in lateral transfer
- Knock Down design is easy to assemble and saves space for transport and storage
- Seat height is adjustable in ½" increments
- Reversible back attachment to accommodate any bathroom
- Attractive platinum grey color
- Maximum Weight Capacity  
350 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



## Specifications:

Overall Width.....	29 ½"
Overall Depth.....	17"
Seat Width.....	25 ½"
Seat Depth .....	16"
Minimum Seat Height.....	19 ¼"
Maximum Seat Height.....	23 ¼"
Clearance Beneath Frame.....	16 ¼" - 20 ¼"
Height of Back Support Above Seat .....	14 ½"
Maximum Weight Capacity .....	350 lb EVENLY DISTRIBUTED



7927KD-1            1 ea  
7927KD-2            2/cs