



- Offers all the features of the Deluxe Forearm Crutches with the added benefit of the therapist-preferred Ortho-Ease[®] handgrip
- Handgrip angle can be adjusted to comfortably fit the user's hand and wrist angle
- The ability to properly position the handgrip can provide correct elbow flexion, which can reduce fatigue and minimize stress on the joints
- Includes laminated color operating instructions
- Attractive bronze color
- Designed to meet the requirements of HCPCS code: E0110
- 300 lb maximum weight capacity - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	<i>Size</i>	<i>Height Adjustment Range</i>	<i>Approximate User Height Range</i>
6345	Large	35" - 41"	5' 10" - 6' 6"
6346	Medium	30" - 36"	5' - 6' 2"
6347	Small	25" - 31"	4' 2" - 5' 2"