

The Lumex Walkabout Basic Four-Wheel Rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 13.5 lbs and supports up to a 250 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.



Features:

- Ergonomic hand grips
- Easy to operate locking and “quick stop” loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable padded backbar with easy release pull knobs
- Adjustable handle height to accommodate broad range of user heights
- Pouch included for storage of personal items
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Three (3) year limited warranty

RJ4900B Blue 1 ea
RJ4900R Red 1 ea



Color options

Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.

Specifications

Width Inside Handgrips.....	16.5"
Overall Width.....	24.2"
Overall Depth.....	24.5"
Seat Width.....	13"
Seat Depth.....	12.5"
Seat Height.....	21"
Caster Diameter.....	6"
Minimum Handle Height.....	32.5"
Maximum Handle Height.....	37.5"
Product Weight.....	13.5 lbs
Weight Capacity.....	250 lb
	EVENLY DISTRIBUTED