

E-Series Hand-Held UV-A Lamps

- Featuring a silver-anodized aluminum housing that makes the lamp both attractive and durable.
- A full range of models offers a selection of UV intensity, size and wattage.
- They are perfect for all applications where a compact lightweight lamp is desired.

1 ea

Model		Typical Peak Intensity ($\mu\text{W}/\text{cm}^2$) at 6"	Dimensions
2211	One 4-Watt integrally filtered tube	800	3 1/4" x 7" x 2 1/2"
2212	One 6-Watt integrally filtered tube	1,280	3 1/4" x 10 1/4" x 2 1/2"
2213	Two 4-Watt integrally filtered tube	620	3 1/4" x 7" x 2 1/2"
2214	One 4-Watt integrally filtered assembly	555	3 1/4" x 8" x 2 1/2"



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com