

**Bariatric Bath Seat**

- Durable blow-molded plastic seat with optional backrest
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- Comes assembled with tool-free back attachment
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0245
- Lifetime limited warranty

*without backrest*

**7932A**                    2/cs

*with backrest*

**7939A**                    2/cs

**7939A-1**                1 ea

<b>Specifications</b> .....	<b>7932A</b>	<b>7939A</b>
Overall Width .....	18"	18"
Overall Depth.....	16 ¾"	20 ½"
Seat Width .....	17 ¾"	17"
Seat Depth.....	12"	16 ¼"
Minimum Seat Height .....	14"	16"
Maximum Seat Height .....	17"	20"
Height of Back Support Above Seat .....	N/A	13"
Maximum Weight Capacity.....	600 lb EVENLY DISTRIBUTED	

