

## Three-In-One Aluminum Commodes

### *With Casters and Removable Back Bar*

- 3-in-1 design is versatile and economical
- Aluminum frame is lightweight and easy to transport
- Contoured arm rests enhance user comfort
- Sturdy plastic seat snaps-on easily
- Comes complete with commode pail, cover and splash shield
- Limited lifetime warranty
- 4" swivel, locking casters
- Designed to meet the requirements of HCPCS code: E0164

**2215B-2**          2/cs

### Specifications

Width Between Arms .....	19 ½"
Overall Width .....	23 ½"
Seat Width .....	13 ¾"
Seat Depth.....	16 ½"
Seat Height.....	21"
Maximum Weight Capacity.....	300 lb EVENLY DISTRIBUTED

