

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame
- Welded steel construction provides strength and stability
- Folds flat to less than 4" for convenient storage and travel
- Tool-free assembly
- Sturdy plastic seat snaps-on easily
- Comes complete with commode pail, cover and splash shield
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0163
- Five (5) year limited warranty



Specifications

| | |
|-------------------------------|---------------------------|
| Width Between Arms | 20" |
| Overall Width | 24" |
| Seat Width | 14" |
| Seat Depth | 16 ½" |
| Seat Height Range | 17" - 23" |
| Maximum Weight Capacity | 300 lb EVENLY DISTRIBUTED |
| Replacement Tip | LX9141 |

7108A 4/cs
7108R Retail Package 2 ea/cs