

Bariatric 3-in-1 Aluminum Commode

- 3-in-1 design is versatile and economical
- Aluminum frame is lightweight and easy to transport
- Contoured arm rests enhance user comfort
- Sturdy plastic seat snaps-on easily
- Comes complete with commode pail, cover and splash shield
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty
- Designed to meet the requirements of HCPCS code: E0168

2190A 2/cs

Specifications

Width Between Arms	22 ¼"
Overall Width	27"
Seat Width	14"
Seat Depth	16 ¾"
Seat Height Range	18" - 23"
Maximum Weight Capacity	400 lb EVENLY DISTRIBUTED
Replacement Tip	LX9141

