

Platinum Collection Three-In-One Steel Commode

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame
- Welded steel construction provides strength and stability
- Increased seat depth provides added stability and comfort
- Durable, white-powder coated steel frame is attractive and easy to clean
- Contoured waterfall armrests are designed to enhance user comfort and provide a good gripping surface
- Sturdy plastic seat snaps-on easily
- Seat height readily adjusts from 17" to 24", in 1 inch increments
- Comes complete with commode pail, cover and splash shield
- Designed to meet the requirements of HCPCS code: E0163
- Limited lifetime warranty

7103 4/cs*Replacement Parts*

- 7103L** Replacement Legs
7103S Replacement Seat and Lid
7103SG Replacement Splash Guard
RP20790-6 Replacement Pail and Lid

Specifications

Width Between Arms	18 ½"
Overall Width	24 ½"
Overall Depth.....	20"
Seat Width.....	14"
Seat Depth.....	16 1/2"
Seat Height Range	17" - 27"
Maximum Weight Capacity.....	300 lb

EVENLY DISTRIBUTED

