

Rolled Foam Mattress

- Economical pressure-reducing and pressure redistribution mattress
- High quality foam mattress offers a high level of therapeutic comfort at an affordable price
- Ships rolled in carton to save on oversize freight charges
- Opens to full size in just minutes
- Zippered, soft, flexible, and breathable Recovery5™ nylon cover moves with the patient to reduce shear and is fluid-resistant, antimicrobial, anti-static and easy to clean
- Meets 16 CFR 1633 requirements
- Designed to meet the requirements of HCPCS code: E0184
- 350 lb maximum weight capacity
- One year limited warranty

GF100 80" x 36" x 6"

