

## Silver Collection 3-in-1 Steel Drop Arm Commode

- Arms drop independently below seat level for safe lateral transfers
- Easy-to-release drop arm mechanism may be operated by patients lacking manual dexterity
- Arms automatically lock in raised position
- Removable backbar facilitates use as a bedside commode, raised toilet seat, or toilet safety frame
- All welded steel construction increases strength and stability
- Durable, grey-coated steel frame is attractive and easy to clean
- Cushioned armrests enhance user comfort
- Sturdy plastic seat snaps-on easily
- Seat height readily adjusts from 18 ¾" to 22 ¾", in 1 inch increments
- Comes complete with commode pail, cover and splash shield
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0165
- Limited lifetime warranty

6433A 1 ea  
 6433A-2 2/cs



**Specifications**

Width Between Arms .....	18"
Overall Width .....	23 ¾"
Seat Width .....	14"
Seat Depth .....	16 ¾"
Seat Height Range .....	18 ¾" - 22 ¾"
Maximum Weight Capacity .....	300 lb EVENLY DISTRIBUTED
Replacement Tip .....	9290