

Home Bed Assist Rail

The Lumex Home Bed Assist Rail offers assistance to the patient/resident for getting in and out of a home style bed.

- Can be used on either side of the bed for assistance with getting in and out of home style bed
- Strong, steel construction with powder-coated finish makes it durable and easy to clean
- Cushioned handle that is adjustable in height assists with getting in and out of bed
- Easily slides in between mattress and bed frame with added feature of anti-slip foam on base bar and securing strap to secure to bed frame
- Extendable legs that touch the floor provide added stability
- Tool-free assembly
- 1 Year Limited Warranty

GF6800 White/Black 1 ea

Specifications

Height adjustment from floor to top of handle	33.5" – 42.5"
Height adjustment from mattress base to top of handle	19" – 21"
Height adjustment from floor to base tubing support under mattress	15" – 22"
Length of rail base tubing support under mattress	43"
Overall width of handle	20"
Product weight	7.5 lb
Weight capacity:	300 lbs

(product is not designed to support the full weight of an individual)

