

Vacutec™ 800 EV2

- Designed for most homecare suction applications
- Low level noise at less than 55 dBA
- Designed for single patient, everyday use
- Easy-to-read gauge and easy-to-turn flow regulator knob
- Suction cup feet for stable operation
- Complies with EN60601 standards
- Complete with disposable 800cc collection bottle
- Automatic shut-off
- Designed to meet the requirements of HPCS code: E0600
- Two years limited warranty

JB0112-016		1 ea
JB0112-016B	Collection Bottle	6/cs
JB0112-016F	Replacement HEPA Filter	3 ea/pk
JB0112-016S	Suction Tubing (one-time patient use) - 6 ft long	
JB0112-016C	Connection Tubing (from Aspirator to Collection Bottle)	

Technical Specifications

Width	4.9"
Depth	15.3"
Height	10"
Weight	11.6 lbs
Power Requirements	115 VAC, 60Hz
Vacuum Range	0 - 560 mm Hg
Power Consumption	60 watts
Operating Flow	0-40 LPM
Bacteria Filter	HEPA
Operating Temperature and Humidity	32°F to 104°F (5° to 40° C)
Operating Humidity.....	95% RH

