

★★★
PATRIOT LX
Full-Electric Homecare Bed



- Robust steel grid deck sleep surface provides superior mattress support and easy to clean surface compared to traditional spring decks
- Head and foot sections can be raised or lowered simultaneously or independently with one hand control
- Recessed main frame offers maximum pinch point protection between grid deck and frame – reducing the risk of injury
- Lightweight design weighs only 126 lbs – lightest bed on the market
- Bed setup is so simple that it can save on setup time and requires minimal maintenance keeping your technicians time efficient
- Patent pending integrated HI/LO mechanism makes bed height adjustment easy and reduces potential loss of parts such as HI/LO rod when used as rentals
- Features dual HI/LO motors to adjust bed height from a low 10" to 20" easy and even allows for trendelenburg and slant positioning
- Patent pending rail slots in frame ensure quick and proper rail placement
- Single motor unit is designed for easy installation, requires less maintenance and provides the ultimate in quiet, smooth head and foot articulation
- Attractive, “home-like” thermo fused bed ends offer a less institutional look and eliminate a bed component that often requires frequent repair and maintenance
- Features SilverSolutions- an anti-microbial technology built into the paint that protects against bacteria and odor
- During a power failure, the motor can lower both the head and foot section with just a single 9 volt battery
- 4" extension kit can be installed on head, foot or both to extend bed length to 84" or 88"

Grid Deck offers the following advantages:

- Much stronger frame
- Better mattress support
- No spring sag or bottoming out with heavier patients
- Easier to clean
- No connecting links or springs to lose

★★★
PATRIOT LX
Full-Electric Homecare Bed



Rail Slot



Casters (optional)



Footboard

Technical Specifications

Total bed weight assembled: 126 lb
 Maximum patient weight: 350 lbs EVENLY DISTRIBUTED
 450 lb distributed total weight capacity
 including patient, mattress, rails and accessories
 Overall Bed Size: 84" L x 35" W
 Lift Range: Head – 60, Knee – 30
 Bed Height: Min. 10", Max. 20" (at deck level)

Warranty:

Bed Frame 5 Years
 Electrical Components 2 Years
 All Other Durable Components 1 Year

<u>Item</u>	<u>Description</u>	<u>HCPCS</u>
US6000	Full-Electric Bed Only	E0266
US6000-IFPKG	Full-Electric Bed with 1633-Innerspring Mattress with Foam Sleep Surface and FDA Full Rails	E0265
US6000-IFPKGHR	Full-Electric Bed with 1633-Innerspring Mattress with Foam Sleep Surface and FDA Half Rails	E0265
US6000-IFPKGQR	Full-Electric Bed with 1633-Innerspring Mattress with Foam Sleep Surface and FDA Quarter Rails	E0265
US6000-PKG	Full-Electric Bed with 1633-Innerspring Mattress and FDA Full Rails	E0265
US6000-PKGHR	Full-Electric Bed with 1633-Innerspring Mattress and FDA Half Rails	E0265
US6000-PKGQR	Full-Electric Bed with 1633-Innerspring Mattress and FDA Quarter Rails	E0265
US6000-PKGC	Full-Electric Bed with 1633-Innerspring Mattress and Chrome Full Rails	E0265
US6000-PKGNR	Full-Electric Bed with 1633-Innerspring Mattress and No Rails	E0265
US6000-RFPKG	Full-Electric Bed with 1633-Rolled Foam Mattress and FDA Full Rails	E0265
US6000-RFPKGHR	Full-Electric Bed with 1633-Rolled Foam Mattress and FDA Half Rails	E0265
US6000-RFPKGQR	Full-Electric Bed with 1633-Rolled Foam Mattress and FDA Quarter Rails	E0265
US6000-RFPKGC	Full-Electric Bed with 1633-Rolled Foam Mattress and Chrome Full Rails	E0265
US6000-RPKG	Full-Electric Bed with No Mattress and FDA Full Rails	E0266
US6000-RPKGHR	Full-Electric Bed with No Mattress and FDA Half Rails	E0266
US6000-RPKGQR	Full-Electric Bed with No Mattress and FDA Quarter Rails	E0266
US6000-SFPKG	Full-Electric Bed with 1633-Standard Foam Mattress & FDA Full Rails	E0265
US6000-SFPKGHR	Full-Electric Bed with 1633-Standard Foam Mattress & FDA Half Rails	E0265
US6000-SFPKGQR	Full-Electric Bed with 1633-Standard Foam Mattress & FDA Quarter Rails	E0265
US6000-SFPKGC	Full-Electric Bed with 1633-Standard Foam Mattress & Chrome Full Rails	E0265