

Neb-u-Lite EV Nebulizer Compressor

- Small footprint and light weight provide exceptional portability
- Compressor weighs only 3.2 lbs
- Reliable, Piston-driven Compressor
- 5 year limited warranty
- Designed to meet the requirements of HCPCS Code: E0570

JB0112-061 12/cs



The John Bunn Neb-u-Lite EV comes complete with:

- 1 Reusable nebulizer (for single patient use only)
- 1 Mouthpiece
- 1 "T" adapter
- 1 PVC tubing
- 3 Spare filters

**Boiling nebulizer components is not recommended*

Technical Specifications

Method Of Nebulization	Piston Driven	Width	5.9"
Power	115VAC, 60Hz	Height	4.5"
Power Consumption	64 watts	Depth	8"
Operating Temperature Range	50°F to 95°F (10°C to 35°C)	Max Pressure (Compressor)	32.5 PSI
Operating Humidity Range	30% to 95% RH	Flow Range (Compressor)	0 - 10 liters/min
Storage Temperature Range	-40°F to 158°F (-40°C to 70°C)	Particle Size (Nebulizer)	0.5 - 5 microns
Storage Humidity Range	10 to 95% RH	Capacity (Nebulizer)	10 ml
Weight	3.6 lb		

** Federal Law (USA) restricts this device to sale by or on the order of a physician. This product should not be used unless the operator has been instructed by a qualified health care professional.*