

Neb-u-Lite® LX2 Nebulizer Compressor

- High-quality, low-noise aerosol compressor with user-friendly design.
- Piston-driven motor ensures high performance and long life.
- Complete with disposable nebulizer featuring 0.2 ml/min flow rate and 5 ml capacity medicine cup and cleaning/instruction manual.
- Designed to meet the requirements of HCPCS code: E0570.
- Five year limited warranty.

JB0112-110 1 ea



Technical Specifications

Method Of Nebulization Piston Driven
 Power 120VAC, 60Hz
 Operating Temperature Range 50°F TO 104°F (10°C TO 40°C)
 Operating Humidity Range 10% to 90% RH non-condensing
 Storage Temperature Range -5°F TO 122°F (-15°C TO 50°C)
 Storage Humidity Range 10% to 90% RH non-condensing
 Weight 4.2 lb (1.91 kg)

Width 6.8"
 Height 4.2"
 Length 11.7"
 Max Pressure (Compressor) 30 psi
 Flow Range (Compressor) 0 - 14 liters/min
 Particle Size (Nebulizer) 0.5 to 5 microns
 Capacity (Nebulizer) 5 ml